ENGLISH

5 Set Range-of-Motion

Depress extension ROM button and slide to desired setting. Depress flexion ROM button and slide to desired setting. Adjust opposite to the same settings of first hinge. **NOTE: Settings automatically lock in place. Refer to standard hinge settings instruction sheet enclosed. CAUTION: The range-of-motion stops on both the medial and lateral hinges must be set the same.**

6 Lock Hinge

Lock hinge at a desired -10°, 0°, 10°, 20°, 30°, and 40° degrees. Bend leg to desired setting and engage lock on each hinge. This will lock the brace in place independent of ROM settings. Make sure ROM settings do not interfere with desired lock setting.

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INDICATIONS: For mild ACL, PCL, and MCL sprains, menisci and patellar retinaculum injuries and mild instabilities.

WARNING: Federal Law Restricts this Device to Sale by or on the Order of a Licensed Health Care Practitioner.

WARNING: For Single Patient Use Only.

CAUTION: Activity modification and proper rehabilitation are an essential part of a safe treatment program. If any additional pain or other symptoms occur while using this device, discontinue use and seek medical attention.

CAUTION: This device is not intended to prevent injury, or reduce or eliminate the risk of re-injury.

CAUTION: The range-of-motion stops on both the medial and lateral hinges must be set the same.

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CAUTION: This device has been developed for the indications listed above, and may be damaged by unusual activities such as sports, falls, or other accidents, or exposure of the hinges to debris or liquids. Do not leave this device in a car on a hot day, or in any other place where the temperature may exceed 150F (65 C.) for any extended period of time. To avoid the risk of accident and possible injury to the patient, discontinue use if such an incident has occurred.

WARRANTY: This device is warranted to be free from defects in material and workmanship for a period of 90 days. These warranties apply to devices that have not been modified or subjected to misuse, abuse or neglect. This device is expressly subject to Bledsoe's Product Line Limited Warranty. Copies of this limited warranty can be obtained from your medical professional, via the internet at www.BledsoeBrace.com, or by calling 1-888-BLEDSOE (253- 3763).

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
CE
Council Directive 93/42/EEC of 14 June 1993 concerning Medical devices


CP020131 Rev J 06/12

© 2012 Medical Technology Inc., d/b/a Bledsoe Brace Systems All Rights Reserved
2601 Pinewood; Grand Prairie, Texas 75051, USA
Toll Free Tel: 1.888.BLEDSOE (253.3763)
Local: 972.647.0884 • Local Fax: 972.660.5495
www.BledsoeBrace.com


EUROPEAN AUTHORIZED REPRESENTATIVE

EC	REP
-----------	------------


Emergo Europe
Molenstraat 15
2513 BH The Hague
The Netherlands
Tel: (+31) 70 345 8570
Fax: (+31) 70 346 7299


2(A)


2(B)


2(C)

2(D)

3

4

ENGLISH

1 Loosen All Straps
Loosen all straps to their maximum circumference.
NOTE: CROSSOVER SHORT has only two straps. *Leave CROSSOVER SHORT D-ring end attached to hinge arm, this allows more freedom when pulling the brace onto leg.*

2 Apply Brace
CLASSIC PULL-ON (A): While seated with the leg bent at 45°, pull the brace onto the leg so that the hinges are lined up with the center of the knee.

WRAPAROUND FC (B): While seated with the leg bent at 45°, wrap the sleeve around the leg so that the hinges are lined up with the center of the knee. Close the lower wrap with the hook-to-pile fastener, re-positioning as needed. Lightly position the hook-to-pile fastener on the upper wrap. If necessary, re-position the wrap for an even fit. The front hole must be positioned directly over the knee cap (patella) as shown **(D)**. **NOTE: Hinge uprights should be checked for alignment prior to securing the straps.**

HYBRID FT (C): While seated with the leg bent at 45°, pull the brace onto the leg so that the hinges are lined up with the center of the knee.

Lightly position the hook-to-pile fastener on the upper wrap. If necessary, re-position the wrap for an even fit. The front hole must be positioned directly over the knee cap (patella) as shown **(D)**. **NOTE: Hinge uprights should be checked for alignment prior to securing the straps.**

CROSSOVER SHORT: While seated with the leg bent at 60°, and without removing the hinges or the straps, pull the sleeve onto the leg until hinges line up with center of knee. The front hole must be positioned directly over the knee cap (patella).

3 Secure Lower Straps
Wrap the lower strap around the brace, then loop the strap through the D-ring. Pull back to tighten, and press to close. Repeat for the remaining lower strap. **NOTE: The CROSSOVER SHORT** has only one lower strap.

4 Secure Upper Straps
Wrap the upper strap around the brace, then loop the strap through the D-ring. Pull back to tighten, and press to close. **NOTE: The CROSSOVER SHORT** has only one upper strap.

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