


Elastic Lumbar Support


102

Warnings and Instructions: Review carefully, proper application is required


⚠ Warning: This device will not prevent or eliminate risk of injury. Do Not Overtighten. If swelling, pain, skin irritation, or an unusual reaction occurs, discontinue use immediately and consult your medical professional.
Care: Hand wash using mild soap. Rinse thoroughly. Air dry only. Do not tumble dry.




1. Center the brace on the spine.



2. Pull and secure front closure.



3. Pull and secure left pull strap. Pull and secure right pull strap.




4. Finished application.

Elastic Lumbar Support


102

Warnings and Instructions: Review carefully, proper application is required


⚠ Warning: This device will not prevent or eliminate risk of injury. Do Not Overtighten. If swelling, pain, skin irritation, or an unusual reaction occurs, discontinue use immediately and consult your medical professional.
Care: Hand wash using mild soap. Rinse thoroughly. Air dry only. Do not tumble dry.




1. Center the brace on the spine.



2. Pull and secure front closure.



3. Pull and secure left pull strap. Pull and secure right pull strap.



4. Finished application.